



Simplicity+		inch	mm
A	Minimum swivel radius	25¾	655
B	Minimum wall to seat back	2½	65
C	Seat depth	16¾	426
D	Minimum folded width to footplate	14½	370
E	Minimum open width to edge of footplate - large footplate (perch seat)	22½ 24¾	570 630
F	Armrest width - external	21¾	550
G	Armrest width - internal	18	460
H	Seat back height	17¾	453
I	Footplate to seat height - standard seat option - perch seat option	17-21¼ 26¾-34½	430-540 677-877
J	Minimum footplate height	2¾	70
Bottom parking space required		22¾ - 27¼ 575mm (50°) to 690 at shallow angle (30°)	
Minimum track intrusion into staircase		7¼	185
Minimum track overhang at top of stairs (typical at 45°)		6	150
Minimum track intrusion at bottom of stairs		9¾	250
Maximum weight limit - hinge/slidetrack		309lb / 140.45 kgs 265lb / 120.45 kgs	
Speed		4¾"/sec (approx)	
Power supply		24v DC battery	
Charger		100-240v (50-60Hz)	
Staircase angle		27-50°	



Sterling Stairlifts from Handicare are backed by a 2 year parts warranty plus a LIFETIME warranty on the gearbox and motor.

Sterling Stairlifts - Handicare Accessibility
 10 Industrial Park Drive
 Mount Pocono, PA 18344
 Toll Free: 866-276-5438
 Local: 570-895-2480
 Fax: 570-895-2481
 info@sterling-stairlifts.com www.sterling-stairlifts.com



Please be aware that actual colors may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.



Simplicity+
 Straight stairlift



The benefits of a Sterling Simplicity+ stairlift from Handicare



If you are thinking about installing a stairlift the Simplicity+ offers you a safe and reliable solution.

Handicare's Sterling Stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Simplicity+ straight stairlift offers you all the basic functions that you would expect plus the option to upgrade to a powered seat swivel and/or folding footplate.

If you have an obstruction at the base of the stairs, you can choose from a manually operated hinged track or the unique SlideTrack that is only available from Handicare.



With a choice of six upholstery colors the Simplicity+ will compliment most color schemes.

The track is fitted to the stairs, not the wall, so you will not need to redecorate following the installation of your Simplicity+ stairlift.

Like all Handicare Sterling stairlifts the Simplicity+ has a battery backup which means it will continue to work in the event of a power outage so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and also comes with two remote controls which means, once folded, you can move it along the track freeing up space in your hallway. The Simplicity+ will charge its batteries anywhere on the track that you choose to leave it.



Safety As Standard

Your Simplicity+ stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

Your Simplicity+ stairlift has an easy release retractable seat belt for added security and peace of mind.

Smooth & Quiet

Because of the high quality motors we use in our Simplicity+ stairlifts we are able to offer perhaps the quietest and smoothest stairlift ride available. All backed by a lifetime gearbox and motor guarantee.

Easy to Use

The operating toggle is easy to use, simply hold it towards the direction you wish to travel, up or down, and your Simplicity+ will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.